

GLOTONE

TOUCH-N-GLOW



FOR MEN & WOMEN

Get all the benefits of Ultrasound, LED, and Thermal therapies in a single, do-it-yourself beauty treatment in the comfort of your home.

All you need to do is just dab your face with water.



HOW DOES GLOTONE WORK?

Glotione combines Ultrasound, LED and Thermal Therapy in beautifying the complexion and at the same time effectively combats skin problems. Each of these therapies individually has been known to provide multiple skin benefits.

Ultrasound

Low Frequency and low intensity ultrasound is used for an anti-bacterial cleansing action that penetrates deep into the skin to increase oxygenation to cells, stimulate collagen protein and aid in the removal of waste products. It benefits the skin by stimulating blood circulation, diminishing fine lines and wrinkles, clarifying acne and blemishes, eliminating irregular pigmentation and redness, and relieving puffiness.

LED Lights

This is a special strong pulse light filling technology, using the colourful wide spectrum light that can rapidly repair and regenerate skin cells, promote collagen, improve skin elasticity, reduce wrinkles, tighten pores, and eliminate pigmentation and acne.

How do the LED Lights of different wavelengths work?

- **Red Light:** (Wavelength 660nm) Brightens and evens the skin tone. Proliferates collagen protein and diminishes the appearance of fine lines and wrinkles. Tightens, firms and repairs damaged skin. An excellent repairing and anti-ageing treatment.
- **Yellow Light:** (Wavelength 590nm) Relieves irritation and aims to enhance immunity within the skin. It also stimulates the lymphatic and nervous system, to rid the skin of toxins and increase circulation. This treatment is best for swollen and sensitive skin.
- **Blue Light:** (Wavelength 475nm) Kills bacterium from within the skin while diminishing inflammation and preventing acne and breakouts.
- **Green Light:** (Wavelength 574nm) Calms and balances the skin tone, and helps to break down and prevent uneven melanin (pigment) production. An excellent treatment that leaves the skin even and bright, while also calming the skin.

Thermal Action

Thermal action deeply cleanses the skin of impurities that can cause a dull and lifeless appearance. Not only does this therapy leave the skin thoroughly clean, but it also improves oxygenation and circulation, and promotes cellular health by staving off free-radical damage and premature ageing. Over time, your skin becomes taught, firm, and bright. Thermal action also helps ultrasound and LED lights to penetrate deeply into the skin, thereby enhancing their effectiveness.

Glotione combines the best aspects of all these therapies in one device to give you a glowing complexion with minimal effort and time.

BENEFITS OF GLOTONE

- **Safe**

Tests have proven Glotone to be 100% safe and highly effective.

1. Conforms to safety standards laid down by IRPA (International Radiation Protection Association) Guidelines Geneva.
2. CE Certified Certificate No. - **TR-LVD- 1609 TCS Turkey European commission.**

- **Can be used on all skin types**

Glotone is safe to be used on all skin types.

- **No known side effects**

No UVA -No UVB-No Ultra-Violet- No Infared.

- **No downtime after treatment**

Glotone is a non-invasive therapy, so you can go about your day normally after use. You can even use makeup immediately after use.

- **Usage**

Operating Glotone is very simple. With no special requirements for a treatment environment, you can use Glotone in the comfort of your home. All you need to do is follow the indicated instructions.

- **Portable**

Glotone comes in a very compact size and is easy to use at home, and carry along when you travel.

- **Low treatment cost**

Compared to photon and laser treatments, Glotone is very cost effective.

Specifications

Model	:	GLOTONE™
Power Supply	:	110VAC/230VAC
Weight	:	995gms

STEPS TO USE GLOTONE

1. Before using the device, wash your face with soap or facewash and dry it.
2. Wet a small piece of cotton with water and dab on your face, or spray some water on your face on the area where you would like to use the device.
3. Connect the cable to the power supply and switch it on.
4. The LED lights on Glotone will start blinking and will also emit a tone indicating that it has started functioning.
5. Hold Glotone by its handle and move it in an upward motion on all areas of your face for one minute.
6. When you are using the device for the first time, the duration can be up to 4 minutes. Subsequently you can use it for 2 minutes or as per your need.
7. You will see the skin starting to glow and the dark circles around the eyes disappearing within a minute.
8. Ladies can use Glotone for a minimum duration of 2 minutes, and since gentlemen have tougher skin, they can use Glotone for a longer duration.
9. You need to use Glotone only twice a week since additional use will not result in any further benefit.

Note: You do not need to use Glotone for more than one session in a day as it will have no further advantage. Glotone will switch off automatically after 5 minutes. And if someone else wants to use Glotone immediately after a session, please restart it. You can clean the diffuser on top of the lights simply with plain water and cotton.

Caution: Never use thinner or other chemicals to clean the diffuser. You may use any disinfectant for cleaning the same.

GLOTONE™
TOUCH-N-GLOW





ISO 9001:2008

ISO 13485:2003



MASER ELECTRONICS (P) LTD.

Corporate Office

Unit No. 307, Lodha Supremus, Road No. 22, Wagle Estate, Near New Passport Office, Thane (W) - 400 604 (India)
Phone +91-22-2158 3100(30 Lines) | Fax : +91-22-2566 5907 | Email : support@glotone.in | Web : www.glotone.in

Registered Office

76, Vishal Industrial Estate, Opp. Ceat Tyres, Bhandup (W), Mumbai - 400078. (INDIA)
Phone +91-22-4111-0000(100 Lines) | Fax : +91-22-2566 5907 | Email : support@maserindia.com | Web : www.maserindia.com